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Use this kit to visualize the concepts of mole and amount of substance. Students will develop a greater understanding of the concept of the mole, as well as a greater knowledge of the physical properties of each element. The kit contains four solids, three bottles for liquids together with 10 balloons for atmospheric air. Common to all substances is that they all consist of 1 mole. The three liquids are not included and have to be supplied at first use. The three liquids are water, ethanol, and acetone, respectively.

## Content

This kit contains 8 different substances for demonstration of the mole concept.

### Four solids

- Carbon
- Aluminum
- Copper
- Bismuth

### Three bottles for liquids

- Water
- Ethanol
- Acetone

### Balloons for gas

- Atmospheric air

## Preparation of the kit


### Liquids

The three bottles are delivered empty and will have to be filled before use. Use a scale to fill the bottles with water, ethanol, and acetone, respectively. Use the table below when filling the bottles.

Substance	Molar mass [g/mol]	Molar volume [mL/mol]
Water	18.0	18.1
Ethanol	46.1	58.4
Acetone	58.1	73.4

## Gas

Inflate a balloon until it contains 24.5 L of air. Assuming the balloon is perfectly round, it will have a circumference of 113 cm.

## Suggestions for extensions of the kit

Substance	Molar mass [g/mol]	Molar volume [L/mol]
Sunflower oil	280.1	0.311
Hydrogen	2.02	22.4
Helium	4.00	24.1
LPG	44-58 <sup>1</sup>	23.3

<sup>1</sup> LPG consists of a mixture of butane and propane. Propane is the lighter with a molar mass of 44 g/mol, whereas butane has a molar mass of 58 g/mol. Regardless of the composition of the mixture the molar volume is 23.3 L/mol.

## Spare parts

151215 Mole balloon, 24 L, 10 pcs.

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© Frederiksen Scientific A/S