Formål med opstillingen

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Purpose of the set-up

This equipment can be successfully used as a demonstration of how vacuum (negative pressure) affects sound and that the air around us has a mass.

The pump plate is used together with the bell jar (178520) to demonstrate the physical condition that occurs for an object / material when a vacuum is formed.

Accessories for the set-up of the pump plate (not included)

- Bell Jar (178510/178520).
- Vacuum pump (069525).
- Wireless Mini Speaker (250520).
- Dasymeter (174010).
- Sound meter (252831).

Technical data

- Diameter, vacuum plates: 206 mm.
- Height: 30 mm.
- Weight: 1.3 kg.

Demostration use

The pump plate is part of experiments where a vacuum is formed.

Two typically used experiments with vacuum take place by examining sound in vacuum. Here, the sound pressure is measured with a

sound meter to demonstrate the change as the vacuum in the bell jar increases.

Another experiment with vacuum is a study of the mass of air. Here it is demonstrated with a dasymeter that when the air is pumped out of the bell jar, the ball on the dasymeter will rise, as the mass of the air disappears.

Installation

The pump plate is connected to a vacuum pump (069525) with a hose via a conical spigot with a diameter of 10 mm. The air in the bell is controlled with a valve tap and a thumbscrew for return air.

In order to create a vacuum, it is required that the compressor is switched on for min. 30 sec..

Experiments

The equipment list for these experiments can be found under the section on "accessories..."

2019-08-14 / MRA

1) **Distribution of sound in vacuum:**

Demonstrate how sound is changed by vacuum.

1. Place a wireless speaker inside the pump bell.
2. Place a sound meter close to the pump bell and switch on the sound meter (dB (A)).
3. Connect PC (mobile phone) via Bluetooth with the speaker and turn on the sound (music).
4. Start the vacuum pump and observe the oscillations on the sound meter.

Demonstrate how the mass of the air in the pump bell is changed by vacuum.

1. Place a dasymeter inside the bell jar and find the balance point horizontally.
2. Start the vacuum pump and observe how the lever on the dasymeter reacts.

In order to create a vacuum, it is required that the compressor is switched on for min. 30 sec..

2) **Air mass in vacuum:**

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© Frederiksen Scientific A/S, Ølgod.