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
This DC power supply provides an output of three different stabilized DC voltages. All outputs are connected to a display where both voltage and current can be read.

The DC outputs are distributed between outputs of 0 to +50 V, 0 to +500 V, and an output of +8 V. All outputs have stepless regulation and they can be regulated independently of each other. They are also galvanically separated by having different connections on the transformer outputs.

Two of the DC outputs can be loaded with up to 50 mA, which requires special attention in use, as this current strength level can be dangerous for exposed persons, e.g. pacemaker. The power supply is also protected against overload by an electric current limiter, and at the same time a warning will be given in the display that the maximum current level has been reached with the text "overload".

### Functionality

Always start by turning the voltage regulator all the way to the left before switching on the mains switch. Repeat this for the other outputs (if used).

If it is experienced that the power supply reports "overload", switch off the power supply at the "red" switch and remove the connected devices.

The power supply can be advantageously used as a voltage source for electro-jet tubes, 508070, and e/m tubes, 507500.

### Specifications

#### 500V

- Output voltage (stabilized): 0-500 V DC
- Infinitely adjustable
- Output current (max.): 50mA
- Ripple and noise (max.): 150 mVpp
- Digital readout: 1% ± 2 LSD

#### 50V

- Output voltage (stabilized): 0-50 V DC
- Infinitely adjustable
- Output current (max.): 50mA
- Ripple and noise (max.): 50 mVpp
- Digital readout: 1% ± 2 LSD

#### 8V

- Output voltage (stabilized): 0-8 V DC
- continuously adjustable
- Output current (max.): 3A
- Ripple and noise (max.): 100 mVpp
- Digital readout: 1% ± 2 LSD

### In general

- Digital readout: 1% ± 2 LSD
- Mains voltage: 230 V AC, 50/60 Hz.
- Dimensions: 312 x 225 x 117 mm.
- Weight: 5 kg

**NB. When using the power supply the following safety rules should be considered:**

#### Caution – HIGH VOLTAGE

- Always work with a partner, never alone
- Wear shoes with rubber soles
- Use only one hand. Keep the other in your pocket or on your back.

This power supply is for teacher use, only.