


# Thermal conductivity of metals



Phy

Physics → Thermodynamics → Temperature & Heat

Physics → Thermodynamics → Heat energy, thermal capacity

Applied Science → Engineering → Materials Science → Thermal & Electrical Properties

Applied Science → Engineering → Renewable Energy → Basic Principles

Applied Science → Engineering → Photonics → Basic Principles



Difficulty level

medium



Group size

2



Preparation time

45+ minutes



Execution time

45+ minutes

**PHYWE**  
excellence in science

# General information

## Application

**PHYWE**  
excellence in science

Fig. 1:  
Experimental  
set-up

Thermal conductivity of metals has many applications in appliances such as heat sinks.

## Other information (1/2)

**PHYWE**  
excellence in science**Prior****knowledge****Main****principle**

The prior knowledge for this experiment is found in the Theory section.

The thermal conductivity of copper and aluminium is determined in a constant temperature gradient from the calorimetrically measured heat flow.

## Other information (2/2)

**PHYWE**  
excellence in science**Learning****objective****Tasks**

The goal of this experiment is to investigate the thermal conductivity of copper and aluminium.

1. Determine the heat capacity of the calorimeter in a mixture experiment as a preliminary test. Measure the calefaction of water at a temperature of 0 °C in a calorimeter due to the action of the ambient temperature as a function of time.
2. To begin with, establish a constant temperature gradient in a metal rod with the use of two heat reservoirs (boiling water and ice water) After removing the pieces of ice, measure the calefaction of the cold water as a function of time and determine the thermal conductivity of the metal rod.

## Safety instructions

**PHYWE**  
excellence in science**Caution:**

Keep the water level such, that the immersion heater is always sufficiently immersed, keep refilling evaporated water during the experiment – the heater will be destroyed by overheating, if the water level is too low.

## Theory (1/2)

**PHYWE**  
excellence in science

If a temperature difference exists between different locations of a body, heat conduction occurs. In this experiment there is a one-dimensional temperature gradient along a rod. The quantity of heat  $dQ$  transported with time  $dt$  is a function of the cross-sectional area  $A$  and the temperature gradient  $dT/dx$  perpendicular to the surface.

$$\frac{dQ}{dt} = -\lambda A \cdot \left( \frac{\partial T}{\partial x} \right) \quad (1)$$

The temperature distribution in a body is generally a function of location and time and is in accordance with the Boltzmann transport equation

$$\frac{\partial T}{\partial t} = \frac{\lambda}{\rho c} \cdot \left( \frac{\partial^2 T}{\partial x^2} \right) \quad (2)$$

Where  $\rho$  is the density and  $c$  is the specific heat capacity of the substance.

## Theory (2/2)


After a time, a steady state

$$\frac{\partial T}{\partial t} = 0 \quad (3)$$

is achieved if the two ends of the metal rod having a length  $l$  are maintained at constant temperatures  $T_1$  and  $T_2$ , respectively, by two heat reservoirs.

Substituting equation (3) in equation (2), the following equation is obtained:

$$T(x) = \frac{T_2 - T_1}{l} \cdot x + T_1 \quad (4)$$



## Equipment

Position	Material	Item No.	Quantity
1	Calorimeter vessel, 500 ml	04401-10	1
2	Calor.vessel w.heat conduct.conn.	04518-10	1
3	Heat conductivity rod, Cu	04518-11	1
4	Heat conductivity rod, Al	04518-12	1
5	Temperature meter digital, 4-2	13618-00	1
6	Universal power supply, 600mA 3/4.5/5/6/7.5/9/12V, incl. 9 adaptors	11078-99	1
7	Immersion probe NiCr-Ni, steel, -50...400 °C	13615-03	1
8	Surface probe NiCr-Ni -50...300°C	13615-04	2
9	Immersion heater, 300W, 220-250VDC/AC	05947-93	1
10	Magnetic stirrer without heating, 3 ltr., 230 V	35761-99	1
11	Magnetic stirring bar 30 mm, cylindrical	46299-02	1
12	Heat conductive paste, 60 g	03747-00	1
13	Gauze bag	04408-00	1
14	Digital stopwatch, 24 h, 1/100 s and 1 s	24025-00	1
15	Tripod base PHYWE	02002-55	1
16	Bench clamp expert	02011-00	1
17	Support rod, stainless steel, 750 mm	02033-00	1
18	Support rod, stainless steel, 1000 mm	02034-00	1
19	Universal clamp	37715-01	4
20	Right angle clamp expert	02054-00	6
21	Supporting block 105x105x57 mm	02073-00	1
22	Beaker, Borosilicate, low-form, 400 ml	46055-00	1
23	Portable Balance, OHAUS CX2200	48921-00	1



# Setup and Procedure

## Procedure (1/5)

1. Measurement of the heat capacity of the lower calorimeter
  - Weigh the calorimeter at room temperature.
  - Measure and record the room temperature and the temperature of the preheated water provided.
  - After filling the calorimeter with hot water, determine the mixing temperature in the calorimeter.
  - Reweigh the calorimeter to determine the mass of the water that it contains.
  - Calculate the heat capacity of the calorimeter.
  - Determine the influence of the heat of the surroundings on the calefaction of the water (0 °C without pieces of ice) by measuring the temperature change in a 30-minute period.

## Procedure (2/5)

**PHYWE**  
excellence in science

### 2. Determination of the thermal conductivity

- Perform the experimental set-up according to Fig. 1.
- Weigh the empty, lower calorimeter.
- Insert the insulated end of the metal rod into the upper calorimeter vessel. To improve the heat transfer, cover the end of the metal rod with heat-conduction paste.
- Attach the metal rod to the support stand in such a manner that the lower calorimeter can be withdrawn from beneath it.
- The height of the lower calorimeter can be changed with the aid of the supporting block. When doing so, care must be taken to ensure that the non-insulated end of the rod remains completely immersed in the cold water during the experiment.

## Procedure (3/5)

**PHYWE**  
excellence in science

- The surface temperature probe must be positioned as close to the rod as possible.
- The outermost indentations on the rod (separation: 31.5 cm) are used to measure the temperature difference in the rod. To improve the heat transfer between the rod and the surface probe, use heat-conduction paste.
- Using an immersion heater, bring the water in the upper calorimeter to a boil, and keep it at this temperature.

**Caution: Keep the water level such, that the immersion heater is always sufficiently immersed, keep refilling evaporated water during the experiment – the heater will be destroyed by overheating, if the water level is too low.**



## Procedure (4/5)

- Ensure that the upper calorimeter is well filled to avoid a drop in temperature due to contingent refilling with water.
- Keep the water in the lower calorimeter at 0°C with the help of ice (in a gauze pouch).
- The measurement can be begun when a constant temperature gradient has become established between the upper and lower surface probes, i.e. when no changes occur during the differential measurement.
- At the onset of measurement, remove the ice from the lower calorimeter.
- Measure and record the change in the differential temperature and the temperature of the water in the lower calorimeter for a period of 5 minutes.
- Weigh the water-filled calorimeter and determine the mass of the water. Settings of the temperature measuring device 4-2:

## Procedure (5/5)

In the first display on the measuring device, the temperature of the lower calorimeter is displayed.

In the second display, the differential measurement between the upper and the lower surface probe is shown.

- The thermal conductivity of different metals can be determined from the measuring results.



# Evaluation

## Task 1

The heat capacity of the calorimeter is obtained from results of the mixing experiment and the following formula:

$$C = c_W \cdot m_W \cdot \left( \frac{\vartheta_W - \vartheta_M}{\vartheta_M - \vartheta_R} \right) \quad (5)$$

$c_W$  = Specific heat capacity of water,  $m_W$  = Mass of the water,  $\vartheta_W$  = Mass of the water,  $\vartheta_W$  = Mixing temperature,  $\vartheta_R$  = Room temperature

The measurement supplies a value of approximately 78 J/K  $\pm$ 25%. The large variations in the results are a result of the manner in which the experiment is performed and of the experimental set-up.

The addition of heat from the surroundings is calculated from the temperature increase  $\Delta T$  of the cold water in the calorimeter.

$$\Delta Q = (c_W \cdot m_W + C) \cdot \Delta T \quad (6) \quad T_0 = \text{Temperature at time } t = 0$$

## Task 2 (1/4)

**PHYWE**  
excellence in science

The heat energy supplied to the lower calorimeter can be calculated using Equation (6). The values and the change in the temperature difference on the metal rod are plotted as a function of time.

In the diagram illustrating the temperature difference, one can see that the temperature essentially remains constant. Consequently, equation 3 can be considered as having been satisfied. In order to calculate the heat energy transported by the metal rod according to Equation 1, the ambient heat fraction must be subtracted.

$$\frac{dQ_{\text{rod}}}{dt} = \frac{dQ_{\text{tot}}}{dt} - \frac{dQ_{\text{sur}}}{dt} \quad (7)$$




Fig. 2: Diagram: Heat of surroundings over time.

## Task 2 (2/4)

**PHYWE**  
excellence in science

$dQ/dt$  for the ambient heat can be calculated from the slope of the graph in Fig. 2.  $dQ/dt$  for the entire set-up can be calculated from the slope of the graph of  $Q$  over  $t$  in Figs. 3 and 4. With the values for the length of the rod ( $\Delta x = 31.5$  cm), the area ( $A = 4.91 \cdot 10^{-4} \text{ m}^2$ ) and the averaged temperature on the metal rod, the heat conduction number can be calculated using Equation (1). The following result as the average values:

$$\lambda_{\text{Al}} = 254 \text{ W/Km}, \lambda_{\text{Cu}} = 447 \text{ W/Km}$$

The literature values are:  $\lambda_{\text{Al}} = 220 \text{ W/Km}$ ,  $\lambda_{\text{Cu}} = 384 \text{ W/Km}$

### Task 2 (3/4)




Fig. 3a:  $\Delta T$  as function of time for aluminium.




Fig. 3b: Q as a function of time for aluminium.

### Task 2 (4/4)




Fig. 4a:  $\Delta T$  as function of time for copper.




Fig. 4b: Q as a function of time for copper.