Koblede harmoniske oscillatorer

27.09.11

Ac 1972.00

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
27.09.11

Ac 1972.00

Coupled harmonic oscillators makes it possible to study harmonic motion on the air track. Up to 5 carts may be linked by springs to study coupled harmonic oscillations or wave propagation.

A single cart can be used with one or more springs to study simple harmonic motion.

The apparatus consists of

- 3 carts for air track
(in addition to the 2 provided with the air track)
- 6 weights for air track carts
(in addition to the 4 provided with the air track)
- 6 springs
- 12 holders for spring
(slightly different from photo above)
- 1 stk. adjustable end stop

Setting up

Set up the apparatus as shown above.
Start with the air track, and align it horizontally as

described in the manual for the air track. Next the carts are placed on the track and connected by springs. The 4 mm plug of the spring connectors fits in the upper of the two holes in the carts as well as in the fixed and adjustable end stops. Turn on the blower and adjust the flow so the carts flow frictionless on the track.

With a 1975.50, 1975.15 or 1975.20 photo gate and a 2002.50 electronic timer you may determine the oscillation period with a suitable flag mounted on the cart. The timer function should be *Period* (see the manual for the timer).

Simple harmonic motion is easily studied with a spring and a single cart by using a string and a weight to tension the spring. Use the pulley supplied with the air track for this. Remember to include the weight in the oscillating mass.

The forks from the rubber band bumpers (provided with the air track) allow one or more springs to be connected to the cart. This way the effect of different spring constants can be studied.

Further investigations

With a 2500.50 or 2501.50 function generator and a 2185.00 vibration generator the coupled oscillators may be driven by an external force down to a frequency of 0,1 Hz.

If a larger amplitude is needed, a ME-8750 mechanical oscillator can be used in the 0,3 – 3,0 Hz range.

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© A/S Søren Frederiksen, Ølgod