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This Geiger-Müller detector is suited for measuring alpha, beta, and gamma radiation.

The detector is provided with a loudspeaker reproducing the signal from the GM tube as loud cracks. The signal is not going into saturation at high intensities.

Moreover, each event is signalled by two flashing blue LEDs.

The detector comprises a 500 V power supply and a GM tube in a common housing.

### Optional accessories

This device is particularly suitable for qualitative demonstration experiments.

However, it is also possible to use the detector with a counter. It is then necessary to acquire a cable:

512560 Cable, modular - Jack 6.3 mm

This cable connects directly to e.g. 513600 GM Counter or Pasco PS-2159 Digital Adapter.

### Specifications

#### Window:

Thickness:	1.5 to 2.0 mg/cm <sup>2</sup>
Effective diameter:	9.14 mm
Material:	Mica

Gas: Neon and halogen

Connection: 6-pin modular socket  
(Used with cable 512560.)

### CAUTION

The GM sensor is provided with a thin and fragile window. The sensor is therefore fitted with a protective cap. The cap is provided with an air vent which must not be covered during mounting and dismantling of the cap. The cap must be removed when measuring alpha and beta radiation, but can be left on when measuring gamma radiation.

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© Frederiksen Scientific A/S