Positiv feltretning

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
Positive field direction

The teslameter is used for measuring magnetic fields up to 2 T. The instrument consists of a detachable probe and a readout box.

The apparatus is powered by a 12 V adapter.

A 10 mm mounting rod for the probe is provided.

The magnetic field is registered with a Hall element that is sensible to magnetic fields *perpendicular to the flat probe tip*.

The probe is provided with an 8-pin DIN plug that also fits Pasco's datalogging equipment (analog adapter PS-2158, or one of the interfaces 500, 550, or 850)

Measures

Probe (Tip):	80 x 8 x 2 mm
Probe (total):	215 x 35 x 25 mm
Readout box:	158 x 108 x 56 mm.

Specifications

Range 100 mT – 1999 mT:	Resolution 1 mT
Range 0.0 – 199.9 mT:	Resolution 0.1 mT

Automatic range selection:

Increasing field: change from 199.9 to 200 mT

Decreasing field: change from 100 to 99.9 mT

Accuracy: 5 %

Calibration when using data logging:

5000 mV ~ 2000 mT

Spare parts

Instruments purchased after ultimo 2011:

355050 Power adapter 12 V, med DC-plug

Instruments purchased before ultimo 2011, use instead :

355035 Power adapter 3-12 V, "universal" (includes Jack plug) – **only** Euro mains plug!

406055 Magnetic field probe

Experiments with the teslameteret

The following descriptions are meant as an inspiration to experiments with magnetic fields. They don't represent complete lab manuals.

Magnetic field from a coil – varying current

Place a coil, e.g. 462520 (400 windings) facing the side of the probe tip. Change the current in 5 to 10 steps, ending at the maximum current allowed for the coil (in the example mentioned, this is 1 A).

Data processing

Plot data in a coordinate system. Describe how the magnetic field depends on the current.

The effect of an iron core

Repeat the above with an iron core in the coil.

Magnetic field from a coil – varying current

Place the probe and a coil so that the coil can be exchanged without changing the geometry.

Select 3 to 4 coils (e.g. 462510, 462520, 462525 and 462530). Find one maximum current allowed for **all** coils and use that for the measurements. Measure the magnetic field from each coil.

Data processing

Plot the magnetic field as a function of the number of windings in a coordinate system. Describe the behaviour.

Magnetic field near the end of a bar magnet

Place a bar magnet (e.g. 330510 – sold in pairs, use only one) on a sheet of millimetre paper to make it easy to find the distance between the end of the magnet and the teslameter probe.

Measure the field along the extension of the magnet at a distance of 10 mm. Repeat in steps of 5 mm until reaching a distance equal to half the length of the magnet.

Plot the results.

Data processing (requires a spreadsheet)

If we consider the bar magnet as consisting of two punctiform magnetic poles, we have measured

the field from one of the poles while the other was "very far away". It will then be natural to guess that the field decreases in proportion to the square of the distance.

Calculate for all distances r the value of $1/r^2$ and plot the magnetic field B as a function of $1/r^2$. Do the data points fall on a straight line?

Try to add a constant distance r_0 to the distances and plot B as a function of $1/(r+r_0)^2$. Vary the size of r_0 until the points lie as close to a straight line as possible. Give an interpretation of the entity r_0 .

Magnetic field from a dipole

Use a neodymium magnet (e.g. 331820 – sold in pairs, use only one).

Use a little adhesive tape to place the magnet on a sheet of millimetre paper to make it easy to find the distance between the end of the magnet and the teslameter probe.

The field is measured along the extension of the centre axis of the magnet.

Start measuring the magnetic field at a distance of 5 mm. Repeat in steps of 5 mm until reaching a total distance of 70 mm. Plot the results.

Data processing (requires a spreadsheet)

As in the experiment above, introduce a correction r_0 . Make a separate column for $r+r_0$.

Plot the measurement data in three log-log coordinate systems: One with all the data, one with data from 5 to 20 mm, one with data from 30 to 70 mm.

In all three cases, fit the data with a power function. Is it reasonable to describe the distance dependency as proportional to $1/r^2$ (with or without r_0)? Other propositions? Comment on the results.

More experiments

The teslameteret is also used in the experiment 138550-EN The beta spectrum - advanced version Find description and lab manual on our web site.

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© Frederiksen Scientific A/S