Intended use

A dip slide culture method for diagnosing urinary tract infections by demonstrating microbes in urine.

Principles of the procedure

The Uricult dip slide system is based on two agar media. One side of the plastic slide is covered with green CLED medium and the other with reddish-brown MacConkey medium for detection of microbes causing urinary tract infections.

The CLED medium is intended for determining the total bacterial count. On the MacConkey medium, bile salts prevent the growth of gram-positive organisms other than enterococci which may grow as pinpoint colonies. This medium supports the growth of gram-negative organisms.

Reagents

Contents

Uricult	Cat. No. 67404
Dip slides	10
Patient labels	10
Instructions for use	1

Typical formulation

CLED medium		MacConkey medium	
Peptone	10.0 g/l	Peptone	20.0 g/l
Meat extract	3.0 g/l	Lactose	10.0 g/l
Lactose	10.0 g/l	Neutral red	0.075 g/l
L-Cystine	0.13 g/l	Bile salts	0.8 g/l
Bromthymol blue	0.03 g/l		

Storage

Uricult is stored at 7...25°C, protected from air and temperature fluctuations. Avoid draughts and storage near heat-generating appliances. **Do not allow to freeze.** The expiry date is marked on the box.

Warnings and precautions

Uricult is for *in vitro* diagnostic use only.

Do not use the product beyond the expiry date marked on the box. Wear protective clothing and disposable gloves while handling samples or tests, and wash hands thoroughly afterwards.

Do not use the Uricult if you detect discoloration or dehydration of the agar, separation of the growth media from the plastic slide or evidence of bacterial or fungal growth.

Because any colonies growing on Uricult are actual or potential pathogens, do not touch the growth.

Sample collection and preparation

Ideally, urine for bacterial culture should remain in the bladder for four hours prior to sampling. Urine samples may be obtained by voiding (clean-voided midstream urine), catheterisation or suprapubic aspiration. The sample should be inoculated onto the Uricult slide immediately after collection. The slide should then at once be returned into its protective tube and the cap closed tightly.

If the urine sample needs to be stored prior to inoculation, it should be maintained refrigerated at 2...8°C no longer than 24 hours.

Uricult test results may be affected if the patient has received anti-infective treatment. The test should not be performed until 48 hours after the final dose of medication.

Procedure

1. Unscrew the slide from the tube without touching the agar surfaces.
2. Holding Uricult by the cap, dip the slide into freshly voided midstream urine so that the agar surfaces are totally immersed. If the volume of urine is too small for this, the agar surfaces can be wetted by pouring urine on them, followed by tilting to ensure complete wetting.
3. Allow excess urine to drain from the slide.
4. Blot the last drops on absorbent paper.
5. Screw the slide tightly back into the tube.
6. Fill in the patient label and attach it to the tube.
7. Place the tube upright in an incubator (36±2°C) for 16–24 hours. The tube may also be sent to a laboratory for incubation.
8. To obtain a colony count (CFU/ml), remove the slide from the tube and compare the colony density with the model chart provided in the kit.

Note:

1. Negative cultures may be incubated for additional 24 hours to ensure that slow-growing bacteria are detected.
2. The inoculated slide may be incubated immediately or stored or transported to a laboratory for incubation and interpretation. Storage or transportation should not exceed 48 hours at 7...25°C, after which Uricult should be incubated at 36±2°C for 16–24 hours. If the slide has been stored or transported for up to 48 hours, only the presence of growth and the colony count should be recorded from it; the colour reaction may be atypical.
3. The inoculated slide may be incubated at room temperature for 1–3 days, after which positive cultures may be sent to a specialised laboratory for further investigation⁵. Negative cultures may be incubated for additional 24 hours to detect slow-growing bacteria⁶.

Quality control

Quality control tests are performed on each lot of Uricult dip slides at the time of manufacture. Should the user wish to perform his own quality control, the following procedure is recommended:

1. Prepare a 10⁵–10⁶ bacteria/ml suspension of each of the following bacteria in sterile saline:
 - a. *Staphylococcus aureus* ATCC 25923
 - b. *Escherichia coli* ATCC 25922
 - c. *Proteus mirabilis* ATCC 12453
2. Use the suspensions to inoculate the Uricult dip slides, using the normal procedure.
3. Interpret the results after a 16–24 hour incubation as follows:

S. aureus ATCC 25923: Growth of colonies on the CLED medium only. Colonies ferment lactose, as indicated by the yellow colour of the colonies and the shift towards yellow of the medium.

E. coli ATCC 25922: Growth of yellow colonies with a shift towards yellow of the CLED medium and growth of pink-red colonies on the MacConkey medium.

P. mirabilis ATCC 12453: Growth of translucent colonies with a shift towards blue of the CLED medium and growth of colourless colonies on the MacConkey medium.

Results' interpretation

After incubation of the inoculated slide, the presence of bacteria is evidenced by colonies on the agar surface. Because a colony is the result of the multiplication of a single bacterial cell, the number of colonies indicates the concentration of colony-forming units (CFUs/ml) in the urine sample. The colony count should be determined from the originally green CLED medium by matching the colony density with the model chart it most closely resembles. It is important to compare the number of colonies, not their size.

The low electrolyte concentration of the CLED medium prevents spreading of *Proteus* strains. Bromthymol blue and lactose in the medium allow the detection of lactose-fermenting bacteria. Such lactose-positive strains grow as yellow colonies and turn the medium yellow, whereas lactose-negative strains grow as translucent colonies with no colour change of the medium.

The originally brownish-red, selective MacConkey medium supports the growth of gram-negative bacteria, but even enterococci may grow as pinpoint colonies on the medium⁷. The selectivity is accomplished by bile salts. Lactose-positive bacteria grow as red and lactose-negative bacteria as translucent colonies on the medium.

When the urinary bacterial content is high (≥ 10⁷ CFU/ml), the agar surfaces may become totally covered by confluent growth. This can be misinterpreted as a negative result. Therefore, any surfaces that appear negative should be examined under a reflecting light; absence of reflection indicates confluent growth. A bright light also allows very small colonies to be detected.

A mixture of different bacterial strains on the Uricult is most likely due to contamination of the urine sample.

Limitations of procedure

Uricult is capable of detecting bacterial concentrations between 10³ and 10⁷ CFU/ml. The model chart allows the determination of colony counts to the nearest power of 10. When the chart is used according to instructions, colony counts show a 99 % correlation with the conventional pour plate method¹.

Expected values

The following values are based on the ECLM-EUG European Urinalysis Guidelines (2000).

Method of sampling clinical status	Significant colony count (CFU/ml)
Midstream, bladder time < 4 hours, symptomatic patient	$\geq 10^3$
Midstream, bladder time > 4 hours	$\geq 10^{4-5}$
Catheter sample from man	$\geq 10^3$
Catheter sample from woman	$\geq 10^4$
Nonsymptomatic bacteriuria	$\geq 10^5$
Puncture sample	Any growth

Note: In some cases, bladder urine < 4 hours may express clinically significant colony counts below 10^3 CFU/ml.


Performance characteristics

Uricult • CLED medium


Arneil, G.C. 1970: Detection of bacteriuria at room temperature. Lancet, January 17, pp 119–121 ⁶		
Number of samples	140	Reference method: Pour plate (Nutrient agar)
Sensitivity	100 %	
Specificity	99 %	
PPV	98 %	
NPV	100 %	

Disposal

Used Uricult dip slides are best disposed of by burning, autoclaving or immersing in a disinfectant overnight, with adherence to local regulations.

Model Chart • Tableau de référence • Standardbildkarte • Tablas de referencia • Tabela de Referência
 Tavola di riferimento • Πρότυπος πίνακας αναφοράς • Modelová tabulka • Referenčna tabela
 Modelkort • Avlesningsmal • Tolkningsmall • Mallitaulu

Literature • Bibliographie • Literatur • Bibliografía • Referências Bibliográficas • Bibliografia • Βιβλιογραφία
Literatura • Referenser • Litteratur • Kirjallisuus

- McAllister TA, Arnell GC, Barr W, Kay P: Assessment of plain dipslide quantitation of bacteriuria. *Nephron* 11: 111–122, 1973.
- Kass EH: Bacteriuria and the diagnosis of infections of the urinary tract. *Archives of Internal Medicine* 100: 709–714, 1957.
- Mackey JP, Sandys GH: Laboratory diagnosis of infections of the urinary tract in general practice by means of a dip-inoculum transport medium. *British Medical Journal* 2: 1286–1288, 1965.
- NCCLS Publication M22-A: Quality Assurance Standards for Commercially Prepared Microbiological Culture Media; Approved standard. Volume 10: 14, 1990.
- Ekman et al.: Streptocult- ja Uricult-levyjen käyttö päivystysaikana. *Aesculapius* 11–12, 1985.
- Arneil GC: Detection of bacteriuria at room temperature. *Lancet*, January 17: 119–121, 1970.
- Granato PA: Evaluation of a dip slide device for enumeration of bacteria in urine. *Laboratory Medicine* Vol. 11, No 4: 246–250, 1980.

Explanation of symbols • Explication des symboles • Zeichenerklärung • Explicación de los símbolos
Explicação de símbolos • Spiegazione dei simboli • Επεξήγηση των συμβόλων • Vysvětlivky použitých symbolů
Pojasnila simbolov • Verklaring van symbolen • Förklaring av symboler • Förklaring av symboler • Symbolien selitykset

IVD
In vitro diagnostic medical device
Dispositif médical de diagnostic *in vitro*
In-vitro-Diagnostikum
Producto sanitario para diagnóstico *in vitro*
Dispositivo médico para diagnóstico *in vitro*
Dispositivo medico-diagnostico *in vitro*
ιατροτεχνολογικό βοήθημα που χρησιμοποιείται για διάγνωση *in vitro*
Diagnostický zdravotnický prostředek *in vitro*
Diagnostická zdravotnícka pomôcka *in vitro*
Medicinsk udstyr til *in vitro*-diagnostik
Medisinsk utstyr for *in vitro*-diagnostikk
Medicinteknisk produkt avsedd för *in vitro*-diagnostik
In vitro -diagnostiikkaan tarkoitettu lääkinällinen laite

REF
Catalogue number
Référence
Bestell Nr.
Número de catálogo
Número de catalogo (n°)
Codice
Αριθμός καταλόγου
Katalogové číslo
Kataloška številka
Bestellingsnummer
Bestillingsnummer
Artikelnummer
Tuotenumero

LOT
Batch code
Número de lot
Lot. Nr.
Número de lote
Número de lote
Numero di lotto
Αριθμός παρτίδας
Číslo šarže
Serijska številka
Batchkode
Lotnummer
Batch nr.
Eräkkoodi

Temperature limitation
Limites de température
Temperaturbegrenzung
Limitación de temperatura
Limites de temperatura
Limiti di temperatura
Περιορισμοί θερμοκρασίας
Teplotní omezení
Teplotné rozmedzie
Temperaturbegrænsning
Temperaturbegrensning
Temperaturbegrensning
Lämpötilarajat

Use by
a utiliser jusqu'à
Verwendbar bis
Fecha de caducidad
Utilizzare até
Utilizzare entro/Scadenza
Χρησιμοποιήστε έως
Použitelné do
Uporabno do
Udløbsdato
Brukes innen
Används före
Käyttävänä ennen

Manufacturer
Fabricant
Hersteller
Fabricante
Produtor
Produttore
Κατασκευάζεται από
Výrobce
Proizvajalec
Fabrikant
Produsent
Tillverkare
Valmistaja

Consult instructions for use
Consulter la notice d'utilisation
Gebrauchsanweisung beachten
Consultense las instrucciones de uso
Consultar as instruções de utilização
Consultare le istruzioni per l'uso
Συμβουλευτείτε τις οδηγίες χρήσης
Viz návod k použití
Sledi navodilu za uporabo
Se bruksanvisningen
Se bruksanvisningen
Läs bruksanvisningen
Katso käyttöohjeita

Sufficient for
Suffisant pour
Ausreichend für
Válido para
Suficiente para
Sufficiente per
Επαρκές για
Dostačuje pro
Zadovoljivo za
Tilstrækkelig til
Tilstrækkelig for
Räcker till
Lukumäärä

Protect from draught and temperature fluctuations
Conserver à l'abri des courants d'air et des fluctuations de température
Vor Zug und Temperaturschwankungen geschützt lagern
Proteger de las corrientes de aire y cambios de temperatura
Proteger das correntes de ar e variações de temperatura
Proteggere da correnti d'aria e variazioni di temperatura
Προστατεύετε από ζήρασια και μεγάλες μεταβολές θερμοκρασίας
Chraňte, před vysušením a kolísáním teploty
Zavarujte pred prepihom in temperaturnimi nihanj
Beskyttes mod tørke og temperatursvingninger
Beskyttes mot trekk og temperatursvingninger
Undvik drag och temperaturvariationer
Suojattava vedolta ja lämpötilan vaihteluilta

Cap • Bouchon • Deckel •
Tapón • Tampa • Tappo •
Καπάκι • Víčko • Pokrov •
Låg • Lock • Kansi

Slide • Lame • Träger •
Lámina • Placa • Piastra •
Πλάκα • Destička • Testna
ploščica • Slide • Levy

Tube • Röhrchen • Tubo •
Cilindro • Σωληναριο •
Nádobka • Epruvetka •
Rør • Rör • Putki

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Uricult® is a registered trademark of Orion Diagnostica Oy.